

UPPER EXTREMITY MEASUREMENT CHART

Bio-Concepts Custom Pressure Garments are available only under Physician's Order



2424 East University Drive, Phoenix, Arizona, U.S.A. 85034-6911
 VOICE: 800-421-5647 / 602-267-7854 • FAX: 800-650-9424 / 602-273-6931
 e-mail: bio-con@bio-con.com • worldwideweb: http://www.bio-con.com

Patient _____
 Last Name First Name
 Date _____

Garment color _____
 Insert _____
 Lining _____

Left arm zipper: Dorsal Palmar Ulnar Radial
 With hook & loop zipper stop tab
 Length (if not full) _____

Right arm zipper: Dorsal Palmar Ulnar Radial
 With hook & loop zipper stop tab
 Length (if not full) _____

Expansion panel: Left Right
 (Expansion panels are full length by default; ulnar when sleeve separate from torso, palmar when attached to a torso garment)

Refer to the Bio-Concepts Measuring Manual for detailed procedures, additional instructions, and example measuring charts.

Additional instructions or comments:

CIRCLE ITEMS ORDERED

Bio-Concepts Use Only: F = _____

Shoulder Flaps

10A. Shoulder Flap regular male & child
 10A. Shoulder Flap regular female
 10B. Shoulder Flap strap male
 10C. Shoulder Flap strap female
 10D. Shoulder Flap strap child

Measure A (armhole circumference)
 with arm down and place the measuring tape distal (away) from the acromion

Male/Child Shoulder Flap Measurements
 Adult Female Shoulder Flap Measurements

	LEFT	RIGHT
A		
B		
C		

Shoulder flaps may be attached to any sleeve that reaches the axilla (No. 1, No. 2, No. 5, No. 8). Select flap style above, sleeve style to lower left. Enter shoulder flap measurements at left and sleeve measurements below.

Sleeves

Sleeves for Gloves & Gauntlets

Sleeves for Briefs, Suits, & Vests

AXILLA
 1 1/2" increments
 WRIST
 Measure with arm down

FROM WRIST or AXILLA	LEFT ARM	RIGHT ARM
0		
1 1/2		
3		
4 1/2		
6		
7 1/2		
9		
10 1/2		
12		
13 1/2		
15		
16 1/2		
18		
19 1/2		
21		

Stump Measurement

STUMP	LEFT	RIGHT
b		
c		

Stump cover style: Orange peel
 Insert